



QF5k Features & Benefits

- Increased max. flow rate of 6,000 lph with most drives
- Self-draining design to minimize nonrecoverable product
- Enhanced venting to reduce the minimum flow rate required to remove entrapped air during priming
- Patented valve plate design to achieve selfdraining and venting
- Up to 120:1 turndown ratio
- Improved linear flow performance
- High flow stability across the entire flow range

- Available in several drive versions:
 - o AC version: With 3-phase asynchronous motor
 - **o Compact drive:** Minimal footprint due to "pump next to motor" design
 - **o HT drive:** Plug-and-play version with integrated motor controller and keypad
 - O-Control: Integrated pump controller with direct sensor connection
- Diaphragm monitoring available as an option
- Motor flange design to reduce pump noise and simplify coupling alignment
- Typical applications include chromatography, TFF, virus filtration, sterile filtration and depth filtration

QF5k Single-Use Quaternary Diaphragm Pumps

- Disposable gamma-irradiated pump chamber
- Pump chamber made of solid polypropylene
- Eliminates the risk of cross-batch contamination
- Lowers cleaning and sterilization costs



QF5k Multiple-Use Quaternary Diaphragm Pumps

- Clean-in-place/steam-in-place (CIP/SIP) and autoclavability
- Stainless steel pump chamber

Quaternary Diaphragm Pumps

Technical Data QF5k

Description			Unit	QF5k-SU	QF5kHT-SU	QF5kQCON-SU	QF5kCD-SU	QF5k-MU	QF5kHT-MU	QF5kQCON-MU	QF5kCD-MU
Flow Rate	Eccentric Shaft 5°	Min.	lph	200 50				200	50		
Flow Rate	Eccentric Shaft 5°	Max.	lph	6,000 5,000			5,000		6,000 5,000		
Pumping Pressure	Media Temperature <40°C	Max.	bar	4				6			
Pumping Pressure	Media Temperature <40°C	Continuous	bar	4				4			
Pumping Pressure	Media Temperature >40°C	Max.	bar	4				6			
Pumping Pressure	Media Temperature >40°C	Continuous	bar	4				4			
Pumping Pressure	Neg. Inlet Pressure	Continuous	bar	-0,2				-0,2			
Temperature		Min.	°C	4°C (39°F)				4°C (39°F)			
Temperature		Max.	°C	60°C (140°F)				60°C (140°F)			
Pump Speed		Min.	rpm	30 13				30		13	
Pump Speed		Max.	rpm	1,200 1,050				1,200 1,050			
Displacement Volume	Eccentric Shaft 5°		ml/rev	91				91			
Connection Specification	Size (inch)		п	1,5"				TriClamp)			
Dimensions	Pump Unit	Length	mm	852	790	893	320	872	851	950	261
Dimensions	Pump Unit	Width	mm	256	275	275	320	257	281	281	320
Dimensions	Pump Unit	Height	mm	333	385	404.5	440	333	385	405	410
Pump Weight with Motor and Housing	Pump Unit		kg	85	100	105	60	95	110	115	70
IP-Protection Class (Total Pump)	Pump Unit			55	54	54	55	55	54	54	55
Filling Volume	Without Connection		ml	788							
Surface	Product Wetted		cm ²	1,586				1,586			
Media	Temperature	Max.	°C	60°C (140°F)				90° (194°F)			
Sterilization	Autoclave	Max.	°C	130°C (266°F)				130°C (266°F)			
Sterilization	Gamma Irradiation	Max.	kGy	50				n/a			
Material	Product Wetted	Pump Housing		PP				1.4435			
Material	Product Wetted	Valve Plate		PP				1.4435			
Material	Product Wetted	Diaphragm		TPE				TPE			
Material	Product Wetted	Valves		EPDM				EPDM			
Material	Product Wetted	O-rings		EPDM				EPDM			



PSG* reserves the right to modify the information and illustrations contained in this document without prior notice. This is a non-contractual document. 11-2022

Data from QF5k-HT or QF5k-QCON drive respectively

QTF-10500-C-01-A4

Authorized PSG® Partner:

Copyright 2022 PSG®, a Dover company

Where Innovation Flows



QUATTROFLOW* Fluid Systems PSG Germany GmbH

PSG Germany GmbH Hochstraße 150-152 47228 Duisburg, Germany Tel: +49 (2065) 89205-0 • Fax: +49 (2065) 89205-40 psgdover.com/quattroflow